

Antibodies Product Datasheet

Product Name:	AWAT2 Antibody
Catalog Number:	C42709
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 260-340
Concentration:	1 mg/ml
Gene Name:	AWAT2 DC4 DGAT2L4 MFAT WS
Protein Name:	Acyl-CoA wax alcohol acyltransferase 2 (EC 2.3.1.75) (Acyl-CoA retinol O-fatty-acyltransferase) (ARAT) (Retinol O-fatty-acyltransferase) (EC 2.3.1.76) (Diacylglycerol O-acyltransferase 2-like protein 4) (Diacylglycerol O-acyltransferase candidate 4) (hDC4
Other Name:	
Human Swiss Prot No.:	Q6E213
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	AWAT2 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year