

Antibodies Product Datasheet

Product Name:	B2LA1 Antibody
Catalog Number:	C42715
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	BCL2A1 BCL2L5 BFL1 GRS HBPA1
Protein Name:	Bcl-2-related protein A1 (Bcl-2-like protein 5) (Bcl2-L-5) (Hemopoietic-specific early response protein) (Protein BFL-1) (Protein GRS)
Other Name:	
Human Swiss Prot No.:	Q16548
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	B2LA1 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year