

Antibodies Product Datasheet

Product Name:	B3A2 Antibody
Catalog Number:	C42717
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 830-910
Concentration:	1 mg/ml
Gene Name:	SLC4A2 AE2 EPB3L1 HKB3 MPB3L
Protein Name:	Anion exchange protein 2 (AE 2) (Anion exchanger 2) (Non-erythroid band 3-like protein) (BND3L) (Solute carrier family 4 member 2)
Other Name:	
Human Swiss Prot No.:	P04920
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	B3A2 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year