

Antibodies Product Datasheet

Product Name:	B3A3 Antibody
Catalog Number:	C42718
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 420-500
Concentration:	1 mg/ml
Gene Name:	SLC4A3 AE3
Protein Name:	Anion exchange protein 3 (AE 3) (Anion exchanger 3) (CAE3/BAE3) (Cardiac/brain band 3-like protein) (Neuronal band 3-like protein) (Solute carrier family 4 member 3)
Other Name:	
Human Swiss Prot No.:	P48751
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	B3A3 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year