

Antibodies Product Datasheet

Product Name:	B3GN1 Antibody
Catalog Number:	C42721
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human B3GN1
Concentration:	1 mg/ml
Gene Name:	B3GNT1 B3GNT6
Protein Name:	N-acetyllactosaminide beta-1,3-N-acetylglucosaminyltransferase (EC 2.4.1.149) (I-beta-1,3-N-acetylglucosaminyltransferase) (iGnT) (Poly-N-acetyllactosamine extension enzyme) (UDP-GlcNAc:betaGal beta-1,3-N-acetylglucosaminyltransferase 1)
Other Name:	
Human Swiss Prot No.:	O43505
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human B3GN1
Storage and Stability:	-20°C/1 year