

Antibodies Product Datasheet

Product Name:	B3GN8 Antibody
Catalog Number:	C42723
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1: 500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human B3GN8
Concentration:	1 mg/ml
Gene Name:	B3GNT8 B3GALT7 BGALT15
Protein Name:	B3GN8
Other Name:	
Human Swiss Prot No.:	Q7Z7M8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of B3GN8 at Human/Mouse
Storage and Stability:	-20°C/1 year