

Antibodies Product Datasheet

Product Name:	B3GT6 Antibody
Catalog Number:	C42725
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 50-130
Concentration:	1 mg/ml
Gene Name:	B3GALT6
Protein Name:	Beta-1,3-galactosyltransferase 6 (Beta-1,3-GalTase 6) (Beta3Gal-T6) (Beta3GalT6) (EC 2.4.1.134) (GAG GalTII) (Galactosyltransferase II) (Galactosylxylosylprotein 3-beta-galactosyltransferase) (UDP-Gal:betaGal beta 1,3-galactosyltransferase polypeptide 6)
Other Name:	
Human Swiss Prot No.:	Q96L58
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	B3GT6 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year