

Antibodies Product Datasheet

Product Name:	B4GT4 Antibody
Catalog Number:	C42728
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	B4GALT4 UNQ552/PRO1109
Protein Name:	Beta-1,4-galactosyltransferase 4 (Beta-1,4-GalTase 4) (Beta4Gal-T4) (b4Gal-T4) (EC 2.4.1.-) (UDP-Gal:beta-GlcNAc beta-1,4-galactosyltransferase 4) (UDP-galactose:beta-N-acetylglucosamine beta-1,4-galactosyltransferase 4) [Includes: N-acetyllactosamine syn
Other Name:	
Human Swiss Prot No.:	O60513
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	B4GT4 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year