

Antibodies Product Datasheet

Product Name:	B4GT7 Antibody
Catalog Number:	C42730
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	B4GALT7 XGALT1 UNQ748/PRO1478
Protein Name:	Beta-1,4-galactosyltransferase 7 (Beta-1,4-GalTase 7) (Beta4Gal-T7) (b4Gal-T7) (EC 2.4.1.-) (UDP-Gal:beta-GlcNAc beta-1,4-galactosyltransferase 7) (UDP-galactose:beta-N-acetylglucosamine beta-1,4-galactosyltransferase 7) [Includes: Xylosylprotein 4-beta-g
Other Name:	
Human Swiss Prot No.:	Q9UBV7
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	B4GT7 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year