

Antibodies Product Datasheet

Product Name:	BACE2 Antibody
Catalog Number:	C42732
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 210-290
Concentration:	1 mg/ml
Gene Name:	BACE2 AEPLC ALP56 ASP21 CDA13 UNQ418/PRO852
Protein Name:	Beta-secretase 2 (EC 3.4.23.45) (Aspartic-like protease 56 kDa) (Aspartyl protease 1) (ASP1) (Asp 1) (Beta-site amyloid precursor protein cleaving enzyme 2) (Beta-site APP cleaving enzyme 2) (Down region aspartic protease) (DRAP) (Memapsin-1) (Membrane-as
Other Name:	
Human Swiss Prot No.:	Q9Y5Z0
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	BACE2 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year