

Antibodies Product Datasheet

Product Name:	ALKB3 Antibody
Catalog Number:	C42744
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 110-190
Concentration:	1 mg/ml
Gene Name:	ALKBH3 ABH3 DEPC1
Protein Name:	Alpha-ketoglutarate-dependent dioxygenase alkB homolog 3 (EC 1.14.11.-) (Alkylated DNA repair protein alkB homolog 3) (DEPC-1) (Prostate cancer antigen 1)
Other Name:	
Human Swiss Prot No.:	Q96Q83
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	ALKB3 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year