

Antibodies Product Datasheet

Product Name:	BC11B Antibody
Catalog Number:	C42751
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 570-650
Concentration:	1 mg/ml
Gene Name:	BCL11B CTIP2 RIT1
Protein Name:	B-cell lymphoma/leukemia 11B (BCL-11B) (B-cell CLL/lymphoma 11B) (COUP-TF-interacting protein 2) (Radiation-induced tumor suppressor gene 1 protein) (hRit1)
Other Name:	
Human Swiss Prot No.:	Q9C0K0
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	BC11B Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year