

Antibodies Product Datasheet

Product Name:	BD1L1 Antibody
Catalog Number:	C42756
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section)
Dilutions:	IHC-p 1:50-300
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 1280-1360
Concentration:	1 mg/ml
Gene Name:	BOD1L1 BOD1L FAM44A KIAA1327
Protein Name:	Biorientation of chromosomes in cell division protein 1-like 1 (Protein FAM44A)
Other Name:	
Human Swiss Prot No.:	Q8NFC6
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	BD1L1 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year