

Antibodies Product Datasheet

Product Name:	BEX3 Antibody
Catalog Number:	C42759
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 1-80
Concentration:	1 mg/ml
Gene Name:	NGFRAP1 BEX3 DXS6984E NADE
Protein Name:	Protein BEX3 (Brain-expressed X-linked protein 3) (Nerve growth factor receptor-associated protein 1) (Ovarian granulosa cell 13.0 kDa protein HGR74) (p75NTR-associated cell death executor)
Other Name:	
Human Swiss Prot No.:	Q00994
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	BEX3 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year