

Antibodies Product Datasheet

Product Name:	BGAT Antibody
Catalog Number:	C42760
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human BGAT
Concentration:	1 mg/ml
Gene Name:	ABO
Protein Name:	BGAT
Other Name:	Histo-blood group ABO system transferase (Fucosylglycoprotein 3-alpha-galactosyltransferase) (Fucosylglycoprotein alpha-N-acetylgalactosaminyltransferase) (Glycoprotein-fucosylgalactoside alpha-N-acetylgalactosaminyltransferase) (EC 2.4.1.40) (Glycoprotei
Human Swiss Prot No.:	P16442
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.47% sodium azide.
Specificity:	This antibody detects endogenous levels of BGAT at Human/Mouse
Storage and Stability:	-20°C/1 year