

Antibodies Product Datasheet

Product Name:	BNIP3L Antibody
Catalog Number:	C42775
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	BNIP3L BNIP3A BNIP3H NIX
Protein Name:	BCL2/adenovirus E1B 19 kDa protein-interacting protein 3-like (Adenovirus E1B19K-binding protein B5) (BCL2/adenovirus E1B 19 kDa protein-interacting protein 3A) (NIP3-like protein X) (NIP3L)
Other Name:	
Human Swiss Prot No.:	O60238
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	BNIP3L Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year