

Antibodies Product Datasheet

Product Name:	BOK Antibody
Catalog Number:	C42779
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 10-90
Concentration:	1 mg/ml
Gene Name:	BOK BCL2L9
Protein Name:	Bcl-2-related ovarian killer protein (hBOK) (Bcl-2-like protein 9) (Bcl2-L-9)
Other Name:	
Human Swiss Prot No.:	Q9UMX3
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	BOK Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year