

Antibodies Product Datasheet

Product Name:	C11B2 Antibody
Catalog Number:	C42803
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 280-360
Concentration:	1 mg/ml
Gene Name:	CYP11B2
Protein Name:	Cytochrome P450 11B2, mitochondrial (Aldosterone synthase) (ALDOS) (EC 1.14.15.4) (EC 1.14.15.5) (Aldosterone-synthesizing enzyme) (CYPXIB2) (Cytochrome P-450Aldo) (Cytochrome P-450C18) (Steroid 18-hydroxylase)
Other Name:	
Human Swiss Prot No.:	P19099
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	C11B2 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year