

Antibodies Product Datasheet

Product Name:	C1QT5 Antibody
Catalog Number:	C42811
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	C1QTNF5 CTRP5 UNQ303/PRO344
Protein Name:	Complement C1q tumor necrosis factor-related protein 5
Other Name:	
Human Swiss Prot No.:	Q9BXJ0
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	C1QT5 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year