

Antibodies Product Datasheet

Product Name:	C2D4D Antibody
Catalog Number:	C42817
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human C2D4D
Concentration:	1 mg/ml
Gene Name:	C2CD4D FAM148D
Protein Name:	C2D4D
Other Name:	C2 calcium-dependent domain-containing protein 4D
Human Swiss Prot No.:	B7Z1M9
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.29% sodium azide.
Specificity:	This antibody detects endogenous levels of C2D4D at Human
Storage and Stability:	-20°C/1 year