

Antibodies Product Datasheet

Product Name:	CA216 Antibody
Catalog Number:	C42828
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human CA216
Concentration:	1 mg/ml
Gene Name:	C1orf216
Protein Name:	UPF0500 protein C1orf216
Other Name:	
Human Swiss Prot No.:	Q8TAB5
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human CA216
Storage and Stability:	-20°C/1 year