

Antibodies Product Datasheet

Product Name:	CAB45 Antibody
Catalog Number:	C42830
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human CAB45
Concentration:	1 mg/ml
Gene Name:	SDF4 CAB45 PSEC0034
Protein Name:	45 kDa calcium-binding protein (Cab45) (Stromal cell-derived factor 4) (SDF-4)
Other Name:	
Human Swiss Prot No.:	Q9BRK5
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human CAB45
Storage and Stability:	-20°C/1 year