

Antibodies Product Datasheet

Product Name:	CAC1G Antibody
Catalog Number:	C42836
Reactivity:	H, Rat
Applications:	IHC-p, IF(paraffin section)
Dilutions:	IHC-p 1:50-300
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 360-440
Concentration:	1 mg/ml
Gene Name:	CACNA1G KIAA1123
Protein Name:	Voltage-dependent T-type calcium channel subunit alpha-1G (Cav3.1c) (NBR13) (Voltage-gated calcium channel subunit alpha Cav3.1)
Other Name:	
Human Swiss Prot No.:	O43497
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	CAC1G Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year