

Antibodies Product Datasheet

Product Name:	CAND1 Antibody
Catalog Number:	C42861
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human CAND1
Concentration:	1 mg/ml
Gene Name:	CAND1 KIAA0829 TIP120 TIP120A
Protein Name:	Cullin-associated NEDD8-dissociated protein 1 (Cullin-associated and neddylation-dissociated protein 1) (TBP-interacting protein of 120 kDa A) (TBP-interacting protein 120A) (p120 CAND1)
Other Name:	
Human Swiss Prot No.:	Q86VP6
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human CAND1
Storage and Stability:	-20°C/1 year