

Antibodies Product Datasheet

Product Name:	CAND2 Antibody
Catalog Number:	C42862
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human CAND2
Concentration:	1 mg/ml
Gene Name:	CAND2 KIAA0667 TIP120B
Protein Name:	Cullin-associated NEDD8-dissociated protein 2 (Cullin-associated and neddylation-dissociated protein 2) (Epididymis tissue protein Li 169) (TBP-interacting protein of 120 kDa B) (TBP-interacting protein 120B) (p120 CAND2)
Other Name:	
Human Swiss Prot No.:	O75155
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human CAND2
Storage and Stability:	-20°C/1 year