

Antibodies Product Datasheet

Product Name:	CAPS2 Antibody
Catalog Number:	C42867
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	CADPS2 CAPS2 KIAA1591
Protein Name:	Calcium-dependent secretion activator 2 (Calcium-dependent activator protein for secretion 2) (CAPS-2)
Other Name:	
Human Swiss Prot No.:	Q86UW7
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	CAPS2 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year