

Antibodies Product Datasheet

Product Name:	CC037 Antibody
Catalog Number:	C42904
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human CC037
Concentration:	1 mg/ml
Gene Name:	C3orf37 DC12
Protein Name:	CC037
Other Name:	UPF0361 protein C3orf37
Human Swiss Prot No.:	Q96FZ2
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.48% sodium azide.
Specificity:	This antibody detects endogenous levels of CC037 at Human/Mouse/Rat
Storage and Stability:	-20°C/1 year