

Antibodies Product Datasheet

Product Name:	CC146 Antibody
Catalog Number:	C42912
Reactivity:	Human, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human CC146
Concentration:	1 mg/ml
Gene Name:	CCDC146 KIAA1505
Protein Name:	CC146
Other Name:	Coiled-coil domain-containing protein 146
Human Swiss Prot No.:	Q8IYE0
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.233% sodium azide.
Specificity:	This antibody detects endogenous levels of CC146 at Human/Rat
Storage and Stability:	-20°C/1 year