

Antibodies Product Datasheet

Product Name:	CCAR1 Antibody
Catalog Number:	C42925
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 840-920
Concentration:	1 mg/ml
Gene Name:	CCAR1 CARP1 DIS
Protein Name:	Cell division cycle and apoptosis regulator protein 1 (Cell cycle and apoptosis regulatory protein 1) (CARP-1) (Death inducer with SAP domain)
Other Name:	
Human Swiss Prot No.:	Q8IX12
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	CCAR1 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year