

Antibodies Product Datasheet

Product Name:	CCD14 Antibody
Catalog Number:	C42928
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human CCD14
Concentration:	1 mg/ml
Gene Name:	CCDC14
Protein Name:	CCD14
Other Name:	Coiled-coil domain-containing protein 14
Human Swiss Prot No.:	Q49A88
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.223% sodium azide.
Specificity:	This antibody detects endogenous levels of CCD14 at Human/Mouse
Storage and Stability:	-20°C/1 year