

Antibodies Product Datasheet

Product Name:	CCG8 Antibody
Catalog Number:	C42946
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 140-220
Concentration:	1 mg/ml
Gene Name:	CACNG8 CACNG6
Protein Name:	Voltage-dependent calcium channel gamma-8 subunit (Neuronal voltage-gated calcium channel gamma-8 subunit) (Transmembrane AMPAR regulatory protein gamma-8) (TARP gamma-8)
Other Name:	
Human Swiss Prot No.:	Q8WXS5
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	CCG8 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year