

Antibodies Product Datasheet

Product Name:	CCL13 Antibody
Catalog Number:	C42948
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 20-100
Concentration:	1 mg/ml
Gene Name:	CCL13 MCP4 NCC1 SCYA13
Protein Name:	C-C motif chemokine 13 (CK-beta-10) (Monocyte chemoattractant protein 4) (Monocyte chemotactic protein 4) (MCP-4) (NCC-1) (Small-inducible cytokine A13) [Cleaved into: C-C motif chemokine 13, long chain; C-C motif chemokine 13, medium chain; C-C motif che
Other Name:	
Human Swiss Prot No.:	Q99616
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	CCL13 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year