

Antibodies Product Datasheet

Product Name:	CCL23 Antibody
Catalog Number:	C42950
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 40-120
Concentration:	1 mg/ml
Gene Name:	CCL23 MIP3 MPIF1 SCYA23
Protein Name:	C-C motif chemokine 23 (CK-beta-8) (CKB-8) (Macrophage inflammatory protein 3) (MIP-3) (Myeloid progenitor inhibitory factor 1) (MPIF-1) (Small-inducible cytokine A23) [Cleaved into: CCL23(19-99); CCL23(22-99); CCL23(27-99); CCL23(30-99)]
Other Name:	
Human Swiss Prot No.:	P55773
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	CCL23 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year