

Antibodies Product Datasheet

Product Name:	CCL25 Antibody
Catalog Number:	C42951
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 30-110
Concentration:	1 mg/ml
Gene Name:	CCL25 SCYA25 TECK
Protein Name:	C-C motif chemokine 25 (Chemokine TECK) (Small-inducible cytokine A25) (Thymus-expressed chemokine)
Other Name:	
Human Swiss Prot No.:	O15444
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	CCL25 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year