

Antibodies Product Datasheet

Product Name:	CCR8 Antibody
Catalog Number:	C42963
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 250-330
Concentration:	1 mg/ml
Gene Name:	CCR8 CKRL1 CMKBR8 CMKBRL2
Protein Name:	C-C chemokine receptor type 8 (C-C CKR-8) (CC-CKR-8) (CCR-8) (CC chemokine receptor CHEMR1) (CMKBRL2) (Chemokine receptor-like 1) (CKR-L1) (GPR-CY6) (GPCY6) (TER1) (CD antigen CDw198)
Other Name:	
Human Swiss Prot No.:	P51685
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	CCR8 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year