

Antibodies Product Datasheet

Product Name:	CD123 Antibody
Catalog Number:	C42969
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 270-350
Concentration:	1 mg/ml
Gene Name:	CDC123 C10orf7 D123
Protein Name:	Cell division cycle protein 123 homolog (Protein D123) (HT-1080) (PZ32)
Other Name:	
Human Swiss Prot No.:	O75794
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	CD123 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year