

Antibodies Product Datasheet

Product Name:	CDK12 Antibody
Catalog Number:	C42981
Reactivity:	Human, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 50-130
Concentration:	1 mg/ml
Gene Name:	CDK12 CRK7 CRKRS KIAA0904
Protein Name:	Cyclin-dependent kinase 12 (EC 2.7.11.22) (EC 2.7.11.23) (Cdc2-related kinase, arginine/serine-rich) (CrkRS) (Cell division cycle 2-related protein kinase 7) (CDC2-related protein kinase 7) (Cell division protein kinase 12) (hCDK12)
Other Name:	
Human Swiss Prot No.:	Q9NYV4
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	CDK12 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year