

Antibodies Product Datasheet

Product Name:	CDK19 Antibody
Catalog Number:	C42986
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 30-110
Concentration:	1 mg/ml
Gene Name:	CDK19 CDC2L6 CDK11 KIAA1028
Protein Name:	Cyclin-dependent kinase 19 (EC 2.7.11.22) (CDC2-related protein kinase 6) (Cell division cycle 2-like protein kinase 6) (Cell division protein kinase 19) (Cyclin-dependent kinase 11) (Death-preventing kinase)
Other Name:	
Human Swiss Prot No.:	Q9BWU1
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	CDK19 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year