

Antibodies Product Datasheet

Product Name:	CE060 Antibody
Catalog Number:	C43005
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human CE060
Concentration:	1 mg/ml
Gene Name:	C5orf60
Protein Name:	CE060
Other Name:	Putative uncharacterized protein C5orf60
Human Swiss Prot No.:	A6NFR6
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.31% sodium azide.
Specificity:	This antibody detects endogenous levels of CE060 at Human
Storage and Stability:	-20°C/1 year