

Antibodies Product Datasheet

Product Name:	CECR1 Antibody
Catalog Number:	C43016
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human CECR1
Concentration:	1 mg/ml
Gene Name:	CECR1 ADA2 ADGF IDGFL
Protein Name:	Adenosine deaminase CECR1 (EC 3.5.4.4) (Cat eye syndrome critical region protein 1)
Other Name:	
Human Swiss Prot No.:	Q9NZK5
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human CECR1
Storage and Stability:	-20°C/1 year