

Antibodies Product Datasheet

Product Name:	CEGT Antibody
Catalog Number:	C43017
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	UGCG
Protein Name:	Ceramide glucosyltransferase (EC 2.4.1.80) (GLCT-1) (Glucosylceramide synthase) (GCS) (UDP-glucose ceramide glucosyltransferase) (UDP-glucose:N-acylsphingosine D-glucosyltransferase)
Other Name:	
Human Swiss Prot No.:	Q16739
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	CEGT Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year