

Antibodies Product Datasheet

Product Name:	CEL3A Antibody
Catalog Number:	C43018
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human CEL3A
Concentration:	1 mg/ml
Gene Name:	CELA3A ELA3 ELA3A
Protein Name:	Chymotrypsin-like elastase family member 3A (EC 3.4.21.70) (Elastase IIIA) (Elastase-3A) (Protease E)
Other Name:	
Human Swiss Prot No.:	P09093
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human CEL3A
Storage and Stability:	-20°C/1 year