

Antibodies Product Datasheet

Product Name:	CENPI Antibody
Catalog Number:	C43022
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 10-90
Concentration:	1 mg/ml
Gene Name:	CENPI FSHPRH1 ICEN19 LRPR1
Protein Name:	Centromere protein I (CENP-I) (FSH primary response protein 1) (Follicle-stimulating hormone primary response protein) (Interphase centromere complex protein 19) (Leucine-rich primary response protein 1)
Other Name:	
Human Swiss Prot No.:	Q92674
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	CENPI Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year