

Antibodies Product Datasheet

Product Name:	CENPK Antibody
Catalog Number:	C43024
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 50-130
Concentration:	1 mg/ml
Gene Name:	CENPK ICEN37 FKSG14
Protein Name:	Centromere protein K (CENP-K) (Interphase centromere complex protein 37) (Protein AF-5alpha) (p33)
Other Name:	
Human Swiss Prot No.:	Q9BS16
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	CENPK Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year