

Antibodies Product Datasheet

Product Name:	CENPR Antibody
Catalog Number:	C43031
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 10-90
Concentration:	1 mg/ml
Gene Name:	ITGB3BP CENPR NRIF3
Protein Name:	Centromere protein R (CENP-R) (Beta-3-endonexin) (Integrin beta-3-binding protein) (Nuclear receptor-interacting factor 3)
Other Name:	
Human Swiss Prot No.:	Q13352
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	CENPR Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year