

Antibodies Product Datasheet

Product Name:	CENPS Antibody
Catalog Number:	C43032
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 70-150
Concentration:	1 mg/ml
Gene Name:	APITD1 CENPS FAAP16 MHF1
Protein Name:	Centromere protein S (CENP-S) (Apoptosis-inducing TAF9-like domain-containing protein 1) (FANCM-interacting histone fold protein 1) (Fanconi anemia-associated polypeptide of 16 kDa)
Other Name:	
Human Swiss Prot No.:	Q8N2Z9
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	CENPS Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year