

Antibodies Product Datasheet

Product Name:	CEP85 Antibody
Catalog Number:	C43037
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human CEP85
Concentration:	1 mg/ml
Gene Name:	CEP85 CCDC21
Protein Name:	CEP85
Other Name:	Centrosomal protein of 85 kDa (Coiled-coil domain-containing protein 21)
Human Swiss Prot No.:	Q6P2H3
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.295% sodium azide.
Specificity:	This antibody detects endogenous levels of CEP85 at Human/Mouse
Storage and Stability:	-20°C/1 year