

Antibodies Product Datasheet

Product Name:	CF062 Antibody
Catalog Number:	C43041
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1: 500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human CF062
Concentration:	1 mg/ml
Gene Name:	C6orf62 XTP12 Nbla00237
Protein Name:	CF062
Other Name:	
Human Swiss Prot No.:	Q9GZU0
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of CF062 at Human/Mouse
Storage and Stability:	-20°C/1 year