

Antibodies Product Datasheet

Product Name:	CH25H Antibody
Catalog Number:	C43057
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human CH25H
Concentration:	1 mg/ml
Gene Name:	CH25H
Protein Name:	Cholesterol 25-hydroxylase (EC 1.14.99.38) (Cholesterol 25-monooxygenase) (h25OH)
Other Name:	
Human Swiss Prot No.:	O95992
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human CH25H
Storage and Stability:	-20°C/1 year