

Antibodies Product Datasheet

Product Name:	CHADL Antibody
Catalog Number:	C43059
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human CHADL
Concentration:	1 mg/ml
Gene Name:	CHADL SLRR4B
Protein Name:	CHADL
Other Name:	Chondroadherin-like protein
Human Swiss Prot No.:	Q6NUI6
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.296% sodium azide.
Specificity:	This antibody detects endogenous levels of CHADL at Human
Storage and Stability:	-20°C/1 year